

## 2X PCR Master Mix

Catalog No.	50501S	50501L
Reactions*	100 reactions	500 reactions

\*Based on 50 µl of reaction volume

### Description

**2X PCR Master Mix** is a ready-to-use 2X master mix solution containing Taq DNA polymerase, dNTPs, MgCl<sub>2</sub> and stabilizers. The master mix allows the setup of the reactions in less than a minute by just adding primers and DNA templates, saves time and reduces contamination due to the fewer pipetting steps required for reaction set up, particularly when preparing large numbers of reactions.

The mix is optimized for efficient and robust amplification of DNA templates with reproducibility and high yields.

Taq DNA Polymerase is derived from the thermophilic bacterium *Thermus aquaticus*. Taq DNA Polymerase is a highly thermostable DNA polymerase that possesses 5'-3' polymerase activities and 5'-3' exonuclease activities, and lacks 3'-5' exonuclease activity. Taq DNA Polymerase produces 3'-dA-tailed PCR fragments that can be used for TA cloning.

The **2X PCR Master Mix** is a reliable, affordable, and robust DNA polymerase for routine PCR applications including genotyping, DNA cloning, library construction, high throughput assays, and screening. Taq DNA Polymerase can amplify DNA target up to 5 kb. The elongation speed is around 1 kb per minute at 72°C.

### Features

- Convenient, ready-to-use 2X master mix for easy reaction setup
- Low background DNA
- Reliable and robust reactions
- Incorporates modified nucleotides such as biotin-, digoxigenin-, and fluorescence-labeled
- PCR fragments have 3'-dA overhangs

### Applications

- Routine PCR
- TA cloning
- High throughput PCR
- Methylated DNA
- Crude sample PCR

### Component

Catalog No.	50501S	50501L
2X PCR Master Mix	2.5 mL	12.5 mL
PCR Water	2.5 mL	12.5 mL

### Storage Condition

- Store kit at -20°C, stable up to 12 months.

### Reagent & Equipment Needed (not provided in this kit)

- Thermal cycler
- PCR plates or PCR tubes
- PCR plate seal film
- PCR grade water

## Master Mix Composition

Component	Concentration
Taq DNA Polymerase	50 U/ml
PCR Buffer	2X
dNTP	0.4 mM each
MgCl <sub>2</sub>	3 mM

## Protocol

### Reaction Setup

Suggested reaction setup (on ice recommended):

Component	Volume	Final conc.
2X PCR Master Mix	25 µl	
DNA template*	variable	1X
Forward Primer (10 µM)	1-2 µl	0.2-0.4 µM
Reverse Primer (10 µM)	1-2 µl	0.2-0.4 µM
Water	variable	
<b>Total Volume</b>	<b>50 µl</b>	

\* Genomic DNA template between 200 ng and 5 ng can be used; cDNA below 100 ng can be used.

### Reaction Condition

Suggested condition:

Step	Temperature	Time	Cycles
Initial denaturation	94°C	60 seconds	1
Denaturation	94°C	30 seconds	25-35 cycles
Annealing	50-65°C	30 seconds	
Extension	72°C	60 seconds**	
Final extension	72°C	2 minutes	1

\*\* Adjust the extension time based on the amplicon size.

### Quality Control

Kit components passed stringent functional quality test.

### Product Use Limitation




This product is developed and sold for research purposes and *in vitro* use only. Please refer to [biodynami.com](https://biodynami.com) for Material Safety Data Sheet of the product.

### Limited Label License

The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynamy is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact [support@biodynami.com](mailto:support@biodynami.com)

### BioDynamy

-  601 Genome Way, Huntsville, Alabama 35806, USA
-  <https://biodynami.com>
-  [support@biodynami.com](mailto:support@biodynami.com)



V1.0 Mar. 2024