

## Library/DNA Size Selection Kit (>30 kb, Magnetic Beads)

Catalog No.	20112S	20112L
Runs*	24	96

\*Based on 50 µl of sample volume


### Description

The **Library/DNA Size Selection Kits (>30 kb, Magnetic Beads)** were developed for DNA size selection using magnetic beads with a selection of DNA fragments over 30 kb. The kits provide a simple and quick approach for the enrichment of a specific range of DNA fragments. The kit workflow allows single-sided size selection for the specific size cutoffs.

The magnetic beads technology uses paramagnetic particles, also known as SPRI (Solid Phase Reversible Immobilization) beads, to bind DNA reversibly and selectively. DNA fragments can be size-selected and purified by changing the properties of the magnetic beads or SPRI beads.

DNA size selection is a selective capture of DNA fragments of a specific range of size for next-generation sequencing (NGS) library preparations, PCR, ChIP assay, DNA ligations, and other applications. Long-read sequencing uses a large DNA fragment as input, the sequencing data is wasted when DNA fragments are short. Usually, DNA size selection is preferred to remove smaller fragments. The kit is ideal for long-read sequencing with DNA size >30 kb.

Specific ranges of NGS libraries can be selected using magnetic beads with different buffer compositions. After the beads-binding step and rinsing step, the NGS library fragments are eluted in water or an appropriate buffer. The magnetic beads method has great advantages over time-consuming column purification and tedious gel-based purification.



Large DNA size selection with single clean-up.

### Features

- High specificity of size selection: >30 kb
- Fast and simple
  - 20-min protocol
  - No gels required
  - No columns required
  - No centrifugation required
- High recovery of selected DNA fragments
- Consistent performance: rapid size selection with high reproducibility

## Component

Catalog No.	20112S	20112L
G6 Beads	1.2 mL	4.8 mL
B30 Buffer	2.0 mL	8.0 mL
WR Buffer	2.0 mL	8.0 mL
Elution Buffer	1.2 mL	4.8 mL

## Storage Condition

- Store at 4°C, stable up to 12 months.

## Reagent & Equipment Needed (not provided in this reagent)

- Magnetic particle concentrator
- 96-well plates or microfuge tubes
- 80% ethanol (prepare before use)

## Protocol

### Attention:

- Invert or shake the beads bottles thoroughly to resuspend the beads before use.
  - Accurate pipetting is needed for precious size selection.
- 1) Add 50  $\mu$ l of the **G6 Beads** to the wells containing 50  $\mu$ l of samples in a 96-well plate. **Slow pipetting of the viscous beads is needed for precise aliquot.** Mix by pipetting gently and thoroughly with a multichannel pipettor. Incubate for 5 min. **Note:** Small sample volumes tend to generate more variations. If the sample volume is less than 50  $\mu$ l, increase the volume to 50  $\mu$ l by adding nuclease-free water.
  - 2) Load the sample plate on a magnet, incubate for 1 min, and discard the supernatant carefully.
  - 3) Remove the plate from the magnet, add 80  $\mu$ l of **B30 Buffer**, mix by pipetting with a multichannel pipettor for 20 times. Make sure the beads is completely resuspended.
  - 4) Add 80  $\mu$ l of **WR Buffer**, mix by pipetting with a multichannel pipettor for 10 times. Incubate for 3 min.
  - 5) Load the sample plate on a magnet, incubate for 2 min, and discard the supernatant carefully.
  - 6) Add 200  $\mu$ l of **80% ethanol** without disturbing the beads. Incubate for 1 min and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
  - 7) Remove the plate from the magnet and resuspend the beads in at least 20  $\mu$ l of water (pH >6.0), **Elution Buffer** (10 mM Tris-HCl), TE buffer, or low TE buffer to elute DNA from the beads. **Note:** Resuspension of the beads in less than 20  $\mu$ l may reduce the yield. A brief centrifugation step may improve bringing eluates to the bottom of the wells before placing on the magnet.
  - 8) Load the plate on the magnet, incubate for 1 min, and transfer the supernatant (containing the size-selected sample) to a new plate or tubes without disturbing the beads.

### Quality Control

Magnetic beads components passed stringent functional quality test.

### Product Use Limitation


This product is developed and sold for research purposes and *in vitro* use only. Please refer to [biodynamy.com](http://biodynamy.com) for Material Safety Data Sheet of the product.

### Limited Label License


The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynamy is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact [support@biodynamy.com](mailto:support@biodynamy.com)

### BioDynamy

 601 Genome Way, Huntsville, Alabama 35806, USA

 <https://biodynamy.com>

 [support@biodynamy.com](mailto:support@biodynamy.com)



Version 1.0 (Apr. 2026)