


Magnetic Beads (DNA & RNA Purification)

Catalog No.	40051S	40051M	40051L
Volume	5 mL	20 mL	100 mL

Description

Solid Phase Reversible Immobilization (SPRI) beads consist of paramagnetic particles coated with carboxyl groups that reversibly bind DNA. They are used for DNA purification because they are fast, simple and efficient. Our **Magnetic Beads (DNA & RNA Purification)** combines BioDynami's proprietary chemistries with reversible DNA-binding properties of SPRI magnetic beads. The beads are developed for effective nucleic acid purification by removing unwanted components such as salts, dNTPs, enzymes, primers, adapters, and other impurities. The Magnetic beads are RNase free, can be used for applications of DNA, and even work with more sensitive RNA without any additional cost.

Our Magnetic beads are optimized to selectively bind DNA fragments of 100 bp and larger, and RNA fragments of 200 bases and larger. Purified DNA and RNA are suitable for downstream applications requiring high quality DNA and RNA, as the purified fragments are free of contaminants and impurities. The beads can be used for NGS library purification, PCR fragment cleanup, molecular cloning, or even nucleic acid concentration.



The beads can also be used for size selection of DNA fragments ranging from 150 bp to 800 bp by changing the bead-to-sample volume ratio and performing single or double-size selection. The beads are an ideal choice for NGS library preparation. They can be easily integrated into the standard workflow of NGS library preparation since the volume ratio is similar for protocols using Magnetic beads.

Features:

- Effective recovery of DNA and RNA samples
 - DNA fragments greater than 100 base pairs
 - RNA fragments greater than 200 bases
- Removal of unwanted components and impurities
- Fragment size selection for specific applications
 - Consistent single or double-size selection
- Flexibility: compatible with manual and automated processing
- Cost effective alternative to AMPure® XP* with equivalent performance

* AMPure® XP is a trade mark of Beckman Coulter.

Applications:

- NGS Library preparation
- Purification of PCR fragments
- Purification of DNA and RNA fragments
- Molecular Cloning
- Other applications requiring purified DNA and RNA

Component

Catalog No.	40051S	40051M	40051L
Magnetic Beads	5 mL	20 mL	100 mL

Storage Condition

- Store at 4°C, stable up to 12 months.

Reagent & Equipment Needed (not provided in this reagent)

- Magnetic particle concentrator
- 96-well PCR plate or Eppendorf tubes
- 80% ethanol (prepare before use)

Protocol

Note: Bring Magnetic beads to room temperature 30 minutes before use. Vortex or invert the bottle to thoroughly resuspend the beads.

Standard DNA & RNA purification

- 1) If sample volume is less than 50 μ l, increase the volume to 50 μ l by adding nuclease-free water. Transfer 90 μ l of the resuspended beads to the sample plate or tube. If sample volume is greater than 50 μ l, scale the beads volume to 1.8X of the sample volume. Mix by pipetting and incubate for 4 min.
- 2) Load the sample plate on a magnet, incubate for 5 min, and discard the supernatant carefully.
- 3) Add 180 μ l of 80% ethanol if the samples are in a 96-well plate; add 800 μ l of 80% ethanol if the samples are in Eppendorf tubes. Incubate for 30 sec and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
- 4) Air-dry the beads on the magnet for 5 min. DO NOT over-dry intact genomic DNA (as it is hard to resuspend it).
- 5) Remove the plate (or tube) from the magnet and resuspend the beads in water or Tris-HCl (10 mM).
- 6) Load the plate on the magnet, incubate for 1 min, and transfer supernatant (containing sample) to a new tube without disturbing the beads.

NGS library purification

- 1) If library volume is less than 50 μ l, increase the volume to 50 μ l by adding nuclease-free water. Transfer 50 μ l of the resuspended beads to the sample plate or tube. If sample volume is greater than 50 μ l, add an equal volume of resuspended beads to the sample. Mix by pipetting and incubate for 3 min.
- 2) Load the plate on a magnet, incubate for 2 min, and discard the supernatant carefully.
- 3) Add 180 μ l of 80% ethanol, incubate for 30 sec, and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
- 4) Remove the plate from the magnet and resuspend the beads in water or Tris-HCl (10 mM).
- 5) Load the plate on the magnet, incubate for 1 min, and transfer supernatant (containing library) to a new tube without disturbing the beads.

Quality Control

Magnetic beads components passed stringent functional quality test.

Product Use Limitation




This product is developed and sold for research purposes and *in vitro* use only. Please refer to BioDynamy.com for Material Safety Data Sheet of the product.

Limited Label License

The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynamy is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact support@BioDynamy.com

BioDynamy

-  601 Genome Way, Huntsville, Alabama 35806, USA
-  <https://BioDynamy.com>
-  support@BioDynamy.com



Jun. 2022