

NGS DNA Library Prep Customization Kit (illumina platform)

Catalog No. 30001S: 24 reactions

Catalog No. 30001L: 48 reactions

Description

The **NGS DNA Library Prep Customization Kit** was developed for construction of high quality libraries for next generation sequencing (illumina platform). The kit adds 3'-dT-tailed library adaptors to both ends of DNA fragments efficiently. The kit needs double strand DNA fragments (blunt and/or sticky) as input DNA for NGS library construction, and is compatible with DNA fragments generated from both enzymatic methods and mechanical methods.

To facilitate the customization of library prep, the kit made it possible for scientists to use their own adaptors, PCR primers and PCR reagents.



Features

- Flexibility: Customized adaptors, PCR primers and PCR reagents allowed
- Fast
 - Total time: <1 hr
 - Hands-on time: 4 min
- Simple work flow: Less reaction steps
- Less magnetic beads required: Reduced more than 50%
- Guaranteed quality: Higher library conversion efficiency
- Input DNA amount: From 100 ng to 1 ug

Component

| | Cat.# 30001S | Cat.# 30001L |
|-------------------------|--------------|--------------|
| • DL1 Buffer | 96 ul | 192 ul |
| • DL1 Enzyme | 96 ul | 192 ul |
| • DL2A Buffer | 576 ul | 1152 ul |
| • DL2 Enzyme | 48 ul | 96 ul |
| • Sodium Chloride (2 M) | 960 ul | 1920 ul |

Storage Condition

- Store kit at -20°C, stable up to 12 months.

Reagent & Equipment Needed (not provided in this kit)

- 3'-end T-overhang adaptor
- PCR primers and PCR reagent (optional)
- Magnetic particle concentrator
- PCR thermal cycler
- 96-well PCR plate
- 80% ethanol (prepare before use)
- Sodium Chloride (2 M)*
- Magnetic Beads (BioDynamix Cat.# 40051) or equivalent

*Sometimes the tube of Sodium Chloride (2 M) may crack during dry ice shipping. Customers need to prepare Sodium Chloride (2 M) in this case.

Protocol

Step 1: End polishing

- 1) Add the following to one well of a 96-well PCR plate:

| | |
|-------------------|---------------------|
| Sheared DNA | 42 ul (100 ng~1 ug) |
| DL1 Buffer | 4 ul |
| <u>DL1 Enzyme</u> | <u>4 ul</u> |
| Total | 50 ul |

- 2) Mix by pipetting ten times.
- 3) Incubate at 20°C for 15 min, 70°C for 5 min. Proceed immediately to step 2.

Step 2: Adaptor addition

- 1) Add the following to Step 1 reaction mixture. Slow pipetting of the viscous DL2 Buffer is needed for precise aliquot.

| | |
|-------------------|-------------|
| Adaptor (20 uM)* | 4 ul |
| DL2A Buffer | 24 ul |
| <u>DL2 Enzyme</u> | <u>2 ul</u> |
| Total | 30 ul |

* Not included in the kit for customization purpose.

- 2) Mix by pipetting ten times.
- 3) Incubate at 20°C for 20 min.
- 4) Add **Sodium Chloride (2 M) 40 ul** to the reaction mixture. Proceed immediately to step 3.

Step 3: Purification

- 1) Resuspend **Magnetic Beads** and transfer 40 ul to the above reaction mixture, mix by pipetting and incubate for 3 min.
- 2) Load the plate on a magnet, incubate for 2 min, and discard the supernatant carefully. Remove all residual supernatant without disturbing the beads.
- 3) Add 180 ul of 80% ethanol, incubate for 30 sec, and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
- 4) Remove the plate from the magnet, resuspend the beads in 22 ul of water.
- 5) Load the plate on the magnet, incubate for 1 min, and transfer 20 ul supernatant (containing library) to a new tube without disturbing the beads.

Step 4: Library amplification when needed (PCR reagent is not provided)

Quality Control

Kit components passed stringent functional quality test.

Product Use Limitation




This product is developed and sold for research purposes and *in vitro* use only. Please refer to BioDynamy.com for Material Safety Data Sheet of the product.

Limited Label License

The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynamy is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact support@BioDynamy.com

BioDynamy

-  601 Genome Way, Huntsville, Alabama 35806, USA
-  <https://BioDynamy.com>
-  support@BioDynamy.com



Oct. 2022