

PCR-free NGS DNA Library Prep Kit (illumina platform)

Catalog No.	30040S	30040L	30041S	30041L
Index type	non-index	non-index	index	index
Reactions	24 reactions	48 reactions	48 reactions	96 reactions

Description

The **PCR-free NGS DNA Library Prep Kit** was developed for construction of high quality libraries for next generation sequencing (illumina platform) without PCR. The PCR-free library can reduce library prep bias, improve coverage of difficult areas, and enhance the detection of sequence variants. The kit is compatible with DNA fragments generated from both enzymatic methods and mechanical methods. Library multiplexing is possible with dual indexes.



Two index types are available for the kit:

Non-index (Cat.# 30040S and 30040L): Libraries do not have index.

Index (Cat.# 30041S and 30041L): Each library contains one i5 index and one i7 index. Library multiplexing up to 96 samples is possible.

Features

- Total time: 1 hr
- Hands-on time: 5 min
- Easy procedure: Ready-to-use master mix & Less reaction components
- Input DNA amount from 100 ng to 1 ug

Component

Catalog No.	30040S	30040L	30041S	30041L
DL1 Buffer	96 ul	192 ul	192 ul	384 ul
DL1 Enzyme	96 ul	192 ul	192 ul	384 ul
DL2A Buffer	576 ul	1152 ul	1152 ul	2304 ul
DL2 Enzyme	48 ul	96 ul	96 ul	192 ul
Sodium Chloride (2 M)***	960 ul	1920 ul	1920 ul	3840 ul
Adaptor	96 ul*	192 ul*	4 ul** X48	4 ul** X96

Storage Condition

- Store kit at -20°C, stable up to 12 months.

Reagent & Equipment Needed (not provided in this kit)

- Magnetic particle concentrator
- PCR thermal cycler
- 96-well PCR plate
- 80% ethanol (prepare before use)
- Sodium Chloride Solution (2 M)***
- Magnetic Beads (BioDynami Cat.# 40051) or equivalent

*non-indexed adaptor

**indexed adaptor

*** Sometimes the tube of Sodium Chloride (2 M) may crack during dry ice shipping. Customers need to prepare Sodium Chloride (2 M) in this case.

Library and Index Information

Non-Indexed library (Cat.# 30040S and 30040L)

Library adaptor:

5'-AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT-3'

5'-GATCGGAAGAGCACACGTCTGAACTCCAGTCACATCTCGTATGCCGCTTCTGCTTG-3'

Indexed library (Cat.# 30041S and 30041L)

Index locations of the indexed libraries:

5' AATGATACGGCGACCACCGAGATCTACACNNNNNNNNACACTCTTTCCCTACACGACGCTCTTCCGATCT-insert-3' TTACTATGCCGCTGGTGGCTCTAGATGTGNNNNNNNNNTGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA-insert-

-insert-AGATCGGAAGAGCACACGTCTGAACTCCAGTCACNNNNNNNNATCTCGTATGCCGCTTCTGCTTG 3'
-insert-TCTAGCCTTCTCGTGTGCAGACTTGAGGTCAGTGNNNNNNNNTAGAGCATACGCAGAAGACGAAC 5'


Note: i5 index: NNNNNNN (in yellow) is the index sequence, 5' to 3' direction.

i7 index: NNNNNNN (in red) is the index sequence, 5' to 3' direction.

List of indexes can be downloaded from:

<https://www.biodynami.com/documents/BioDynami-PCR-free-Index.xls>

For Cat.# 30041S and 30041L, the indexed adaptors have been aliquot into 96-well plates by column (i5) and row (i7). Each well contains one i5 index and one i7 index as shown in figure. Example: library #71 has i5-11 index and i7-06 index.



Protocol

Step 1: End polishing

- 1) Add the following to one well of a 96-well PCR plate:

Sheared DNA	42 ul (100 ng~1 ug)
DL1 Buffer	4 ul
<u>DL1 Enzyme</u>	<u>4 ul</u>
Total	50 ul

- 2) Mix by pipetting ten times.
- 3) Incubate at 20°C for 15 min, 70°C for 5 min. Proceed immediately to step 2.

Step 2: Adaptor addition

- 1) Add the following to Step 1 reaction mixture. Slow pipetting of the viscous DL2A Buffer is needed for precise aliquot.

DL2A Buffer	24 ul
Adaptor	4 ul
<u>DL2 Enzyme</u>	<u>2 ul</u>
Total	30 ul

- 2) Mix by pipetting ten times.
- 3) Incubate at 20°C for 20 min.
- 4) Add **Sodium Chloride (2 M) 40 ul** to the reaction mixture. Proceed immediately to step 3.

Step 3: Beads purification

- 1) Resuspend **Magnetic Beads** and transfer **40 ul** to the above reaction mixture, mix by pipetting and incubate for 3 min.
- 2) Load the plate on a magnet, incubate for 2 min, and discard the supernatant carefully. Remove all residual supernatant without disturbing the beads.
- 3) Add 180 ul of 80% ethanol, incubate for 30 sec, and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
- 4) Remove the plate from the magnet, resuspend the beads in in **22 ul of water** or **Tris-HCl (10 mM)**.
- 5) Load the plate on the magnet, incubate for 1 min, and transfer 20 ul supernatant (containing library) to a new tube without disturbing the beads.

Quality Control

Kit components passed stringent functional quality test.

Product Use Limitation

This product is developed and sold for research purposes and *in vitro* use only. Please refer to BioDynami.com for Material Safety Data Sheet of the product.

Limited Label License




The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynami is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact support@BioDynami.com

About PCR master mix:

This product is licensed from Bio-Rad Laboratories, Inc. under U.S. Pat. Nos. 6,627,424, 7,541,170, 7,560,260, 7,670,808, 7,666,645, 7,919,296, 8,232,078, 8,367,376, 8,415,129, 8,445,249, 8,470,573, 8,476,045, 8,895,283, and 8,900,846 and corresponding patents in other countries for use only in: (a) standard (non-real time) PCR in the research field only, but not digital PCR; (b) real-time PCR for internal product research and development purposes only, and not for sales to end-users within the research field; (c) any in-vitro diagnostic application, including applications using real-time PCR, but not digital PCR; and (d) any non-PCR applications in DNA sequencing, isothermal amplification, and the production of synthetic DNA.

BioDynami

 601 Genome Way, Huntsville, Alabama 35806, USA
 <https://BioDynami.com>
 support@BioDynami.com



Oct. 2022