

RNA Contamination Removal Magnetic Beads (gDNA Purification)


Catalog No.	40060S	40060L
Runs*	40 runs	200 runs

*Based on sample volume of 100 μ l



Description

Genomic DNA extraction is the most popular laboratory technique. However, it is common that the extracted genomic DNA is contaminated with RNA. To reduce the impact on downstream applications, it is critical to purify the genomic DNA by removing RNA contamination.

We have developed a simple reagent to completely remove RNA contamination in the extracted genomic DNA samples based on Solid Phase Reversible Immobilization (SPRI) beads. The magnetic beads are coated with carboxyl groups that can reversibly bind to DNA in an appropriate buffer. Our **RNA Contamination Removal Magnetic Beads** combines the reversible DNA-binding properties of SPRI magnetic beads with BioDynamix's unique chemistries. The reagent removes the contaminated RNA and recovers the genomic DNA in the same reaction step. Moreover, unwanted components such as proteins, dNTPs, salts, enzymes, and other impurities can also be removed in the same step.



The genomic DNA can be used for downstream applications such as PCR, QPCR, sequencing, microarray, and enzymatic treatment etc. The beads can be an effective and inexpensive reagent for RNA depletion for routine genomic DNA purification.



Features

- Effective removal of RNA contamination by RNase
- High recovery rate of genomic DNA by magnetic beads
- Simple and fast beads-based protocol
- Removal of unwanted components and impurities

Component

Catalog No.	40060S	40060L
RNA Contamination Removal Magnetic Beads	40 runs	200 runs

Storage Condition

- Store at 4°C, stable up to 12 months.

Reagent & Equipment Needed (not provided in this reagent)

- Magnetic particle concentrator
- 96-well PCR plate or 1.5 ml tubes
- Thermal cycler, incubator, or heat block
- 80% ethanol (prepare before use)

Protocol

Note: Invert or shake the bottle to thoroughly resuspend the beads.

- 1) Transfer samples to a 96-well plate or 1.5 ml tubes, and add 0.5X volumes of the beads to the wells containing samples. **Slow pipetting of the viscous beads is needed for precise aliquot.** Mix by pipetting gently and thoroughly.
- 2) Incubate at 37°C for 10 min. **Note:** Incubation time can be extended in case of heavy RNA contamination.
- 3) Load the sample plate on a magnet, incubate for 3 min, and discard the supernatant carefully.
- 4) Add 180 µl of 80% ethanol without disturbing the beads. Incubate for 2 min and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads. Repeat 80% ethanol wash one more time.
- 5) Air-dry the beads on the magnet for 1-2 min.
- 6) Remove the plate from the magnet and resuspend the beads in at least 20 µl of water, Low TE Buffer or TE Buffer. **Note:** Resuspending the beads in less than 20 µl will reduce the yield.
- 7) Load the plate on the magnet, incubate for 1 min, and transfer the supernatant (containing DNA) to a new tube without disturbing the beads.

Quality Control

Magnetic beads components passed stringent functional quality test.

Product Use Limitation

This product is developed and sold for research purposes and *in vitro* use only. Please refer to biodynami.com for Material Safety Data Sheet of the product.

Limited Label License

The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynamy is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact support@biodynami.com

BioDynamy

-  601 Genome Way, Huntsville, Alabama 35806, USA
-  <https://biodynami.com>
-  support@biodynami.com



Aug. 2023