

## miRNA Purification Kit (tRNA Depletion)

Catalog No.	40057S	40057L
Runs*	24 runs	96 runs

\*Based on 50 µl of sample volume


### Description

The **miRNA Purification Kit (tRNA Depletion)** is developed for microRNA (miRNA) purification by removal of tRNA, 5S RNA, and other larger RNA contaminations.

microRNAs are typically around 22 nucleotides long and can be used in various downstream applications relating to gene expression regulation, miRNA profiling, and functional analysis. However, miRNA isolated from commercial microRNA extraction kits usually contain contaminations of tRNA, 5S RNA, and other larger RNA molecules. The contamination may result in decreased sensitivity or inconsistent data for applications such as PCR, RT-PCR, next-generation sequencing (NGS), microarray, etc.

We have developed a Solid Phase Reversible Immobilization (SPRI) beads-based technology to solve the problem. SPRI beads are well used for nucleic acid purification. The paramagnetic beads coated with carboxyl groups can reversibly bind to nucleic acid that are 100 base pairs or longer. However, with traditional SPRI beads, nucleic acids shorter than 100 bp are poorly recovered.

Our proprietary magnetic bead technology not only overcomes the short nucleic acid recovery hurdle, but also enables effective separation of miRNA from large RNA molecules. With this technology, the kit selectively purifies miRNA and removes tRNA, 5S RNA, and other larger RNA molecules. The magnetic beads can also remove impurities and unwanted components such as salts, proteins, dNTPs, and detergents.



### Features

- Purification of miRNA by removing tRNA, 5S RNA, and larger RNA contaminations
- Removal of impurities and unwanted components

### Component

Catalog No.	40057S	40057L
S1 Beads	1.2 mL	4.8 mL
EM Buffer	1.8 mL	7.2 mL
S2 Beads	0.24 mL	0.96 mL

### Storage Condition

- Store at 4°C, stable up to 6 months.

### Reagent & Equipment Needed (not provided in this reagent)

- Magnetic particle concentrator
- 96-well PCR plate or microcentrifuge tubes
- 80% ethanol (prepare before use)

### Protocol

**Note:** Invert or shake the beads to thoroughly resuspend the beads.

- 1) Transfer samples (up to 50  $\mu$ l) to a 96-well plate.
- 2) Add 1X volumes of the **S1 Beads** to the wells containing samples (see Reagent volume calculation table below). **Slow pipetting of the viscous beads is needed for precise aliquot.** Mix by pipetting gently and thoroughly. Incubate for 3 min.
- 3) Load the sample plate on a magnet, incubate for 3 min, and transfer the supernatant to a new plate.
- 4) Add 1.5X volumes of **EM Buffer** to the supernatant and mix by pipetting (see Reagent volume calculation table below).
- 5) Add 10  $\mu$ l of **S2 Beads** and mix by pipetting, incubate 10 min.
- 6) Load the sample plate on a magnet, incubate for 10 min, and discard the supernatant carefully.
- 7) Add 180  $\mu$ l of 80% ethanol without disturbing the beads. Incubate for 1 minute and discard the supernatant carefully. Remove all residual ethanol without disturbing the beads.
- 8) Air-dry the beads on the magnet for 3 min.
- 9) Remove the plate from the magnet and resuspend the beads in water or Tris-HCl (10 mM).  
**Note:** Resuspending the beads in less than 20  $\mu$ l may reduce the yield.
- 10) Load the plate on the magnet, incubate for 1 min, and transfer supernatant (containing microRNA) to a new tube without disturbing the beads.

**Reagent volume calculation table**

Sample volume	S1 Beads volume (1X of sample volume)	EM Buffer volume (1.5X of sample volume)
30 $\mu$ l	30 $\mu$ l	45 $\mu$ l
40 $\mu$ l	40 $\mu$ l	60 $\mu$ l
50 $\mu$ l	50 $\mu$ l	75 $\mu$ l

### Quality Control

Magnetic beads components passed stringent functional quality test.

### Product Use Limitation


This product is developed and sold for research purposes and *in vitro* use only. Please refer to BioDynami.com for Material Safety Data Sheet of the product.

### Limited Label License

The product is developed and sold exclusively for research purposes and *in vitro* use only. The product or its any individual component has not been tested for use in diagnostics or drug development, and is not suitable for administration to human or animal.

The purchaser of this product is granted a limited, non-transferable right to use the purchased amount of the product only for internal, research purposes for the sole benefit of the purchaser. The buyer cannot sell or otherwise transfer (i) this product (ii) its components or (iii) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. This product is for internal research purposes only and is not for use in commercial purposes of any kind. "Commercial purposes" includes any activity for which a party receives consideration and may include, but is not limited to, (1) use of the product or its components or derivatives in manufacturing, (2) transfer or sale of vectors made with the product or components or derivatives of the product, (3) use of this product or components or derivatives of the product made therefrom to provide a service, information, or data to a third party in return for a fee or other consideration, or (4) resale of the product or its components or derivatives, whether or not such product or its components or derivatives are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, BioDynami is willing to accept return of the products with a full refund. For information on obtaining additional rights, please contact [support@biodynami.com](mailto:support@biodynami.com)

## BioDynami

 601 Genome Way, Huntsville, Alabama 35806, USA

 <https://biodynami.com>

 [support@biodynami.com](mailto:support@biodynami.com)



May 2023