

DNA Concentrator (Magnetic Beads)

The Easy & Fast Way to Concentrate DNA/RNA



The DNA Concentrator (Magnetic Beads) can be used to efficiently concentrate DNA/RNA from various samples with low concentration without the need of a DNA vacuum concentrator. Optimized to concentrate DNA and RNA of 20 bp and larger.

Features

- Concentrate DNA and RNA samples, even small DNA, RNA, oligos
- Binding capacity: 10 ug DNA per prep
- Remove unwanted components and other impurities
- No vacuum
- No centrifuge
- No column

Low Cost

\$1.4 / sample

Concentrate



DNA/RNA:
20 bp & larger

No vacuum

No centrifuge

No column

Recovery Efficiency



Product Name	Cat.#
DNA Concentrator (Magnetic Beads)	40055

More details at <https://biodynami.com/product/dna-concentrator/>



Web: biodynami.com
Tel: (256) 872-0805
Email: support@biodynami.com

Fast
Simple
Affordable
High-Quality

