

# DNA Ladder


BioDynami DNA ladders from 20 bp to 20 kb.  
Ready-to-load



## 20 bp DNA Ladder

### Features


- Composed of eight bands as shown below
- Premixed with 6X DNA loading buffer for direct gel loading



## 100 bp DNA Ladder

### Features


- Composed of nine bands as shown below. The 500-bp band with higher concentration is easily distinguishable from the others
- Premixed with 6X DNA loading buffer for direct gel loading



## 1 kb DNA Ladder

### Features


- Composed of ten bands as shown below. The 4 kb band with higher concentration is easily distinguishable from the others
- Premixed with 6X DNA loading buffer for direct gel loading



## 1 kb Plus DNA Ladder

### Features

- Composed of 13 bands as shown below. The 10 kb and 4 kb bands with higher concentration are easily distinguishable from the others
- Premixed with 6X DNA loading buffer for direct gel loading



### Product Name

20 bp DNA Ladder  
100 bp DNA Ladder  
1 kb DNA Ladder  
1 kb Plus DNA Ladder

### Cat.#

10002  
10003  
10004  
10005

More details at: <https://biodynami.com/dna-ladder/>



Web: [biodynami.com](http://biodynami.com)  
Tel: (256) 872-0805  
Email: [support@biodynami.com](mailto:support@biodynami.com)

Fast  
Simple  
Affordable  
High-Quality

