

# Magnetic Beads (Short Oligo Purification)

## The Way to Purify as Short as 6 nt Oligos

The **Magnetic Beads (Short Oligo Purification) Kit** enables the recovery of oligos as short as 6 nt. 80% of the 6 nt oligos and 90% of the >8 nt oligos can be recovered.



The **Only** beads purifying short oligo

Works with oligos

**>6 nt**

### Features

- Effective purification of short oligos
  - 6 nt oligos: 80% recovery rate
  - >8 nt oligos: >90% recovery rate
- Removal of impurities and unwanted reaction components



Product Name	Cat.#
Magnetic Beads (Short Oligo Purification)	40058

More details at <https://biodynami.com/product/short-oligo-purification/>



Web: [biodynami.com](http://biodynami.com)  
Tel: (256) 872-0805  
Email: [support@biodynami.com](mailto:support@biodynami.com)

Fast  
Simple  
Affordable  
High-Quality

