

cfDNA Purification Kit (Magnetic Beads)



Separate cfDNA and gDNA; Recover both



The cfDNA Purification Kit (Magnetic Beads) was developed for cfDNA purification by separating genomic DNA contamination from cfDNA. The kit is effective for recovery of both cfDNA and genomic DNA separately.

Features




- Separate cfDNA from genomic DNA, and recover both types of DNA
- Recover cfDNA effectively
 - As short as 50 bp can be recovered
- Recover high molecular weight genomic DNA effectively
- Remove unwanted components and other impurities



Separate
cfDNA from gDNA

Recover
cfDNA

Recover
gDNA (optional)



Kit list

- cfDNA Purification Kit (Magnetic Beads), 24 runs
- cfDNA Purification Kit (Magnetic Beads), 96 runs

Cat.#

- 40056S
- 40056L

More details at <https://biodynami.com/product/cfdna-purification-magnetic-beads/>



Web: biodynami.com
Tel: (256) 872-0805
Email: support@biodynami.com

Fast
Simple
Affordable
High-Quality

